

Little Anne QCPR



Bystander CPR saves lives

QCPR now available in Little Anne.

Providing feedback based on watching your students perform CPR can be challenging. The next generation of Little Anne includes QCPR measurement and feedback technology that takes the guesswork out of instruction – all in the palm of your hand.

All existing Little Anne manikins can be upgraded with the QCPR technology, enabling:

- Real-time feedback on compressions and ventilations
- A performance summary with an overall score and improvement suggestions
- A detailed overview of compression release, depth and rate, ventilation volume, number of compressions/ventilations/cycles

QCPR technology has been proven to be intuitive and effective for performing and improving quality CPR – helping save more lives.



www.laerdal.com



Laerdal
helping save lives



QCPR Instructor app

WHY

- Sets focus on quality in community CPR training
- Facilitates better community CPR training
- Quality CPR training is strongly associated with improved sudden cardiac arrest survival
- Objective measurement is fun and engaging
- Makes it easy, powerful and affordable to measure and improve community CPR training

HOW

Allows the instructor to:

- quickly identify who needs personal assistance
- provide accurate and tailored instructions
- provide objective assessment and informative feedback
- ensure substantial hands-on practice to meet skill performance objectives

WHAT

- QCPR feedback on individual performance
- Realtime and summative feedback available on tablet and TV/monitor for up to six manikins simultaneously

QCPR Learner app



052017V01 ©2017 Laerdal Medical. All rights reserved.



SkillGuide



QCPR Learner app
QCPR Instructor app



QCPR Learner app



Laerdal
helping save lives

For more information, visit www.laerdal.com